



Stellar M22

Shine on



One modular system. A galaxy of treatments.

4 high-end technologies can treat more than 30 skin conditions, including:

Photorejuvenation, pigmentation, vascular lesions, rosacea, scars, textural issues, and hair removal. As Stellar M22™ is a modular system, you have the flexibility to add technologies as your practice expands.



Stellar™ IPL

A Single handpiece with changeable filters and lightguides.



ResurFX™

The only true fractional non-ablative laser, with CoolScan™ scanner.



Multi-Spot™ Nd:YAG

For effective treatment of leg veins and vascular lesions.



Q-Switched Nd:YAG

A real 6-8 nsec pulse with homogeneous beam profile ensures perfect treatment control for skin toning by treatment of pigmented lesions, and dark tattoos removal.

Before & After Results

Before



4 months after 2 treatments



Photofractional™ treatment (Sequential treatment of IPL + ResurFX™), In-house clinic

Baseline




After 3 treatments



Treatment with IPL Vascular filter. Courtesy of Matteo Tretti Clementoni, MD

The Stellar M22™ OPT® Advantages

- ◆ Control pulse shape by equalizing energy distribution over the entire pulse or sequential pulses
- ◆ Gentler, more comfortable, patient-friendly procedures, with lower effective fluences
- ◆ High peak power, shorter pulses - optimal for IPL skin treatments using photorejuvenation and treatment of benign pigmented lesions
- ◆ Advanced OPT will also allow determining the specific fluence per sub-pulse when using multiple-sequential pulsing (MSP™), for fine-tuned treatment settings



Multiple Sequential Pulsing (MSP™)


- ◆ Protects the epidermis, allowing the safe use of higher fluences
- ◆ Improves the efficacy in clearance of vascular lesions*

How it works:

Light energy is delivered during the pulses, and the tissue cools during the delays. More energy is safely delivered to the target.

What it does:

MSP™ allows for cooling between pulses to safely remove hair on all skin types including darker skin and reduces the chance of skin damage in vascular lesion treatments.



ResurFX™ CoolScan™ Scanning Shapes

6 different shapes, with sizes and densities changeable from the user interface, allow for treatment precision and control.



Intuitive User Interface

With easy to learn treatment presets. The Stellar M22 is easy to both learn and use, with the Lumenis treatment presets, and is ergonomically designed with practitioners in mind.

Quick replacement of handpieces / filters / light guides is enabled for fast and customizable treatments, while large storage compartments, both on the system and in the cart enable a seamless experience.



Warnings and risks: Stellar M22 is contraindicated for patients with active infection, viral, fungal or bacterial diseases. The use of Stellar M22 could cause redness, swelling, scarring, change of pigmentation. See the system's user manual for a complete list of contraindications and risks.

Stellar M22, Stellar IPL, OPT, MSP, ExpertFilters, SapphireCool lightguides are trademarks or registered trademarks of the Lumenis Group of Companies. Specifications are subject to change without notice.

Lumenis Ltd.
Yokneam Industrial Park
6 Hakidma Street P.O.B. 240
Yokneam 2069204, Israel
Tel: +972-4-959-9000

Lumenis (Germany) GmbH
Heinrich-Hertz-Str 3 D-63303
Dreieich-Dreieichenhain
Germany
Tel: +49 (0) 6103 8335 0

EMEA
Dreieich Dreieichenhain,
Germany
T +49 6103 8335 0
F +49 6103 8335 300

Roma, Italy
T +39 06 90 75 230
F +39 06 90 75 269

Hertfordshire, UK
T +44 20 8736 4110
F +44 20 8736 4119

JAPAN
Tokyo, Japan
T +81 3 6743 8300
F +81 3 6743 8301

AMERICAS
San Jose, CA, USA
T +1 408 764 3000
+1 877 586 3647
F +1 408 764 3999

ASIA PACIFIC
Gurgaon, India
T +91 124 422 07 95

Kowloon, Hong Kong
T +852 217 428 00
F +852 272 251 51

Sydney, Australia
T +61 1800 5863 64